DAC-SpectroMed s.r.l.

MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64
Tel.: /+37322/ 574900, 574922/23; факс: /+37322/ 574920
Email: office@dacspectromed.com
www.dacspectromed.com

PT MD 11-38623324-002:2002

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PT MD 11-38623324-002:2002

ALP - DAC


ALKALINE PHOSPHATASE
KINETIC TEST WITH DIETHANOLAMINE BUFFER

For «in vitro» diagnostic use only
Store at 2-8°C

Cod 2005A100 100 ml

PRINCIPLE

Alkaline phosphatase (ALP) catalyzes in alkaline medium the transfer of the phosphate group from 4-nitrophenylphosphate to diethanolamine (DEA), liberating 4-nitrophenol. The catalytic concentration is determined from the rate of 4-nitrophenol formation, measured at 405-410 nm^{1,2}.

CONTENTS AND COMPOSITION

| | | |
|------------------------|-----------|------------|
| Reagent A | 2 x 50 ml | pH 10,2 |
| Diethanolamine | | 1,0 mol/l |
| Magnesium chloride | | 0,6 mmol/l |
| Reagent B | 2x0,225 g | |
| 4-nitrophenylphosphate | | 10 mmol/l |

Harmful! Harmful if swallowed

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

Indications of deterioration:

Reagents: presence of particulate material, turbidity, absorbance of the blank over 1,30 at 405(±10) nm (1 cm cuvette).

SAMPLES

Serum free of Hemolysis.

Alkaline Phosphatase in serum or plasma is stable for 7 days at 2-8°C.

REFERENCE VALUES

Children: 180-1200 U/l.

Adult: 98-279 U/l.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 405-410 nm.

Cuvettes with 1 cm light path.

Pipettes for 20 µl and 1,0 ml.

PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious.

REAGENT PREPARATION

Working Reagent:

0,225 g Reagent B + 50 ml Reagent A. Mix gently.
Working Reagent is stable for 2 weeks at 2-8°C.

PROCEDURE

Assay conditions

| | |
|---------------|----------------------|
| Wavelength: | 405-410 nm |
| Temperature: | 37°C |
| Cuvette: | 1 cm light path |
| Read against: | distilled water |
| Method: | kinetic (increasing) |

1. Bring the Working Reagent and the instrument to reaction temperature (37°C).
2. Pipette into a cuvette 1 cm light path:
Working Reagent 1,0 ml
Sample or Standard 20 µl
3. Mix and insert the cuvette into the photometer. Start the stopwatch.
4. After 60 sec record initial absorbance and record absorbance at 1 minute intervals thereafter for 2 minutes.
5. Calculate the difference between consecutive absorbances, and the average absorbance difference per minute (ΔA/min).

CALCULATIONS

The ALP concentration in the sample is calculated using the following general formula:

$$\frac{\Delta A / \text{min}_{\text{Sam}}}{\Delta A / \text{min}_{\text{St}}} \times C_{\text{St}} = C_{\text{Sam}}$$

Calculation using factor
405 nm: Activity (U/l) = ΔA / min_{Sam} x 3000

METROLOGICAL CHARACTERISTICS

Detection limit: 3,0 U/l.

Linearity limit: 1800 U/l.

Repeatability (within run):

| Mean Concentration | CV | n |
|--------------------|--------|----|
| 191 U/l | 0,91 % | 20 |
| 642 U/l | 0,96 % | 20 |

Reproducibility (run to run):

| Mean Concentration | CV | n |
|--------------------|--------|----|
| 176 U/l | 2,29 % | 25 |
| 416 U/l | 1,43 % | 25 |

CV - coefficient of variation n - number of determinations

Interferences: Bilirubin 128,3 µmol/l (0,075 g/l), lipid 10 g/l, glucose 55,5 mmol/l (10 g/l) and ascorbic acid 2,84 mmol/l (0,5 g/l) do not interfere. Other drugs and substances may interfere⁴.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Alkaline Phosphatase catalyzes the hydrolysis of organic phosphate monoesters at alkaline pH. The enzyme is present in practically all tissues of the body, especially at or in the cell membranes, and it occurs at particularly high concentrations in placenta, intestinal epithelium, kidney tubules, osteoblasts and liver. The form present in the sera of normal adults originates mainly in the liver and bone.

Elevated serum ALP is found in patients with bone disease associated with increased osteoblastic activity (Paget's disease, primary and secondary hyperparathyroidism, bone tumors, rickets, osteomalacia, bone fractures) and also in patients with hepatobiliary disease (obstructive jaundice, hepatitis, hepatotoxicity caused by drugs, liver cancer). Physiological changes, such as bone growth and pregnancy, may cause increases in ALP levels^{5,6}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

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