DAC-SpectroMed s.r.l.

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Cod 3008B500 500 ml

PRINCIPLE

The bilirubin in the sample reacts with diazotized sulfanilic acid with forming of diasobilirubin and colored complex. The intensity of coloration, measured at 500-560 nm or 595 (± 10) nm, is proportional to bilirubin concentration.

The bilirubin contained in serum (plasma) is conventionally divided in two fractions: "direct" (conjugated, water soluble) and "indirect" (uncombined, insoluble). Direct bilirubin determination uses water medium, and total (direct+indirect) uses caffeine reagent as diluent.

CONTENTS AND COMPOSITION

Reagent 1	2x5,0 g (2x200 ml)	
Caffeine		25.75 mmol/l
Sodium benzoate		52 mmol/l
Sodium acetate		92 mmol/l
Reagent 2	0.4 g	
Sulfanilic acid		
Reagent 3	6 ml	
Hydrochloric acid		2.45mol/l
Reagent 4	12.5 mg	
Sodium nitrite		
Reagent 5	9.0 g	
Sodium hydroxide		

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label. Caffeine Reagent, Diazoreagent I and Diazoreagent II are stable for 3 months at 2-8°C. Working Reagent 5 is stable for 3 months at 15-30°C in polyethylene bottle.

SAMPLES

Serum (plasma).

Bilirubin in serum (plasma) is stable within 3 hours at 15-30°C, 1 day at 2-8°C in dark place.

ADDITIONAL EQUIPMENT AND REAGENTS

Analyzer, spectrophotometer or photometer with filter 520-560 nm or 595 (± 10) nm. Thermostat 60°C. Timer. Bilirubin Standard. Droppers for 100 μ l, 250 μ l, 500 μ l and 1,75 ml

PRECAUTION

The kit is for in vitro use only.

All specimens must be considered potentially hazardous and handled as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

REAGENT PREPARATION

Caffeine Reagent: dissolve contents of the bottle with Reagent 1 quantitatively with 150 ml of the distilled water in volumetric flask, dilute with the distilled water warming up to 50-60°C and mix thoroughly. After cooling, bring volume to 200 ml with the distilled water.

Diazoreagent I: dissolve, while heating, the contents of the bottle with Reagent 2 quantitatively with 60 ml of the distilled water in flask. Add the contents of the bottle with Reagent 3 and mix thoroughly. Only after complete dilution of sulphanyllic acid and cooling of the solution, bring the volume to 80 ml with distilled water and mix thoroughly.

Diazoreagent II: dilute contents of Reagent 4 in 2,5 ml of distilled water.

Diazomixture: must be prepared immediately before use. Mix 10 ml of Diazoreagent I with 0,3 ml of Diazoreagent II, or any other volumes based upon ratio 100:3 v/v.

Working Reagent 5: dilute contents of Reagent 5 bottle in 30 ml of distilled water.

PROCEDURE

1. Pipette into labeled test tubes:

Sample, ml	Direct bilirubin		Total bilirubin	
	Blank	Sample	Blank	Sample
0,5 ml	0,5 ml	0,5 ml	0,5 ml	0,5 ml
Caffeine Reagent, ml	-	-	1,75 ml	1,75 ml
Physiological solution	2,0 ml	1,75 ml	0,25 ml	-
Diazomixture, ml	-	0,25 ml	-	0,25 ml
Incubate at 16-25°C	5-10 min	5-10 min	20 min	20 min
Working Reagent 5	0,1 ml	0,1 ml	0,1 ml	0,1 ml

2. Read the absorbance (A_S) of the Sample and Blank (A_B) at 595(± 10) nm against the distilled water.

Coloration is stable for at least 60 minutes.

CALCULATIONS

For the calculation of direct and total bilirubin concentrations in sample, calculate the difference between absorbencies of Sample (A_S) and Blank (A_B) for direct (DB) and total (TB) bilirubin.

Using these values, and calibration curve, plotted under total bilirubin standards, determine values of total and direct bilirubin.

REFERENCE VALUES

Adults²:

Total	Up to 1,1 mg/dl = 18,8 μ mol/l
Direct	Up to 0,25 mg/dl = 4,3 μ mol/l

Newborns³:

Age	Premature	Full-term
Up to 24 h	1,0-8,0 mg/dl = 17-137 μ mol/l	2,0-6,0 mg/dl=34-103 μ mol/l
Up to 48 h	6,0-12 mg/dl=103-205 μ mol/l	6,0-10 mg/dl=103-171 μ mol/l
3-5 days	10-14 mg/dl = 171-239 μ mol/l	4,0-8,0 mg/dl = 68-137 μ mol/l

These values are approximate, it is recommended that each laboratory should establish its own reference values.

QUALITY CONTROL

It is recommended to use the Control Serum level I and II (Sera N-DAC Cod 2055S5 and Sera P-DAC Cod 2057S5) to verify the performance of the measurement procedure.

METROLOGICAL CHARACTERISTICS

Sensitivity limit:	0.03 mg/dl = 0.51 μ mol/l.
Linearity limit:	9.92 mg/dl = 170 μ mol/l.

For higher values dilute sample 1/2 with saline solution and repeat measurement.

Coefficient of variation not more than 7 %.

Interference:

Hemoglobin (10 g/l) does not interfere. Lipemia (triglycerides > 5 g/l) interferes. Other drugs and substances may interfere⁴.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

Bilirubin is the product of disintegration of heme moiety of the hemoglobin released from aged or damaged erythrocytes, which are destroyed in the reticuloendothelial cells. Bilirubin then is transported to the liver in association with albumin. Inside the hepatocytes bilirubin is conjugated with glucuronic acid and is excreted into bile. A number of inherited and acquired diseases affect production, uptake, metabolism, and excretion of bilirubin, resulting in hyperbilirubinemia^{3,5}. Unconjugated hyperbilirubinemia is observed in newborns (physiological jaundice), in increased red cell destruction (hemolytic anemia, extensive hematoma), in ineffective erythropoiesis and in some rare genetic diseases. Conjugated hyperbilirubinemia is associated with decreased excretion of bile due to liver diseases or to intrahepatic or extrahepatic cholestasis.

NOTES

Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard (Multi S1-DAC cod 2051M5).

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