DAC-SpectroMed s.r.l.

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PT MD 11-15796482-001:2003

Crea Kinetic–DAC.Lq

CREATININE
 KINETIC METHOD JAFFE WITH ALKALINE PICRATE

For « in vitro » use only
 Store at 15-30°C

Cod 3035C100	100 ml
Cod 3035C500	500 ml
Cod 3035C1000	1000 ml

PRINCIPLE

Creatinine in the sample reacts with picrate in alkaline medium forming a colored complex. The intensity of coloration, measured at 500(±10) nm, is proportional to creatinine concentration.

The intensity of complex formation is measured in a short period to avoid interferences^{1,2}.

CONTENTS AND COMPOSITION

Reagent A	
Picric acid	25 mmol/l
Reagent B	
Sodium hydroxide	0.2 mol/l
Detergent	
Irritant! Avoid the ingress on skin and mucous membrane.	
Standard	5 ml
creatinine Aqueous primary standard	2 mg/dl (177 µmol/l)

STORAGE AND STABILITY OF REAGENTS

Store reagents at 15-30°C and use until the expiry date shown on the label. In case of precipitation in Reagent B it is recommended to warm at 37±20C till completely dissolved prior to use.

Indications of deterioration: absorbance of the blank over 0,350 at 500 nm (1 cm cuvette).

SAMPLES

Serum (plasma) or urine.
 Dilute fresh urine 1/50 with distilled water before measurement.
 Heparin, EDTA, oxalates and fluorides may be used as anticoagulants.
 Creatinine in samples is stable for 24 hours at 2-8°C.

REFERENCE VALUES

Serum and Plasma: Men: 0.9 – 1.3 mg/dl = 80 - 115 µmol/l
 Women: 0.6 – 1.1 mg/dl = 53 - 97 µmol/l
 Urine: Men: 14 - 26 mg/kg/24-h = 124-230 µmol/l/kg/24-h
 Women: 11-20 mg/kg/24-h=97-177 µmol/l/kg/24-h

These values are approximate, each laboratory should establish its own reference values.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with filter 500(±10) nm.
 Thermostat for 37°C. Dropper for 100 µl and 1.0 ml. Timer.

PRECAUTION

The kit is only for in vitro use. All specimens must be considered potentially hazardous and handled as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

QUALITY CONTROL

It is recommended to use the Control Serum N end P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action when controls do not recover within the acceptable tolerances.

REAGENT PREPARATION

Creatinine Standard is provided ready to use.
 Working Reagent: Mix equal volumes of Reagent A and Reagent B and let stand for 30 minute at room temperature (16-25°C).
 Stable for 7 days at 2-8°C.

PROCEDURE

Method: kinetic (increasing)
 Wavelength: 500 (±10) nm
 Temperature: 37°C
 Read against: Working Reagent

- Bring the Working Reagent and the photometer to reaction temperature (37°C).
- Pipette into a cuvette:

Working Reagent, ml	1.0
Sample or Standard (St), µl	100
- Mix and insert the cuvette into the photometer. Start the timer.
- Record absorbance at 500(±10) nm after 30 seconds (A₁) and after 90 seconds (A₂).
- Calculate the difference between consecutive absorbencies
 Sample – (A₂ - A₁)_{Sample} and Creatinine Standard – (A₂ - A₁)_{St}.

CALCULATIONS

The creatinine concentration (K_{Sample}) in the sample is calculated using the following general formula:

$$\frac{(A_2 - A_1)_{Sample}}{(A_2 - A_1)_{St}} \times K_{St} \times K_F = K_{Sample}$$

If the Creatinine Standard provided has been used to calibrate:

	Serum and Plasma	Urine
$\frac{(A_2 - A_1)_{Sample}}{(A_2 - A_1)_{St}} \times 2 =$ mg/dl creatinine		$\times 100 =$ mg/dl creatinine
$\frac{(A_2 - A_1)_{Sample}}{(A_2 - A_1)_{St}} \times 177 =$ µmol/l creatinine		$\times 8840 =$ µmol/l creatinine

METROLOGICAL CHARACTERISTICS

Sensitivity limit: 0.03 mg/dl = 2.65 µmol/l
 Linearity limit: 20 mg/dl = 1768 µmol/l

At higher concentration dilute sample with distilled water 1:2 and repeat the measurement.

Repetition (within period):

Mean concentration	CV*	n*
1,7 mg/dl = 150 µmol/l	2,9 %	20
5,3 mg/dl = 468 µmol/l	1,3 %	20

Reproducibility: (from period to period):

Mean concentration	CV*	n*
1,7 mg/dl = 150 µmol/l	3,9 %	25
5,3 mg/dl = 468 µmol/l	2,9 %	25

* CV - coefficient of variation n – number of determinations

Sensitivity: 31 ΔmA. dl/mg = 0,351 ΔmA. l/µmol.

Interference: bilirubin (10 mg/dl), hemoglobin (0,1 g/l), protein and ketonic bodies do not interfere. Lipemia (triglycerides > 2 g/l) interfere. Other substances and drugs may interfere⁴.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

Creatinine is a catabolic end product of creatine (phosphocreatine). The amount produced each day is related to the muscle mass. Creatinine is freely filtered by the glomerulus (small amount is reabsorbed and then is secreted by the renal tubules). Creatinine measurement is used exclusively in the assessment of kidney function (impaired renal perfusion, loss of functioning nephrons) and in the monitoring renal dialysis^{3,5}.

Clinical diagnosis must be made on integration of clinical and laboratory data.

BIBLIOGRAPHY

- Bartels H, Bohmer M. Eine mikromethode zur kreatininbestimmung. *Clin Chim Acta* 1971; 32: 81-85.
- Fabiny DL, Ertingshausen G. Automated reaction-rate method for determination of serum creatinine with CentrifChem. *Clin Chem* 1971; 17: 696-700.
- Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press