## DAC-SpectroMed s.r.l.

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PT MD 11-15796482-001:2003

## Fructosamine-DAC

FRUCTOSAMINE  
KINETIC-PHOTOMETRIC

For « in vitro » use only  
Store at 2-8°C

Cod 3037F100 2x50 ml

### PRINCIPLE

Serum glycated proteins reduce nitroblue tetrazolium (NBT) salts in alkaline medium.

The intensity of the color measured at 530(±20) nm is proportional to the glycated proteins in the sample.

### CONTENTS AND COMPOSITION

Reagent 2x50 ml pH 10,35  
NBT 0,25 mmol/l  
Carbonate buffer 0,1 mmol/l

Fructosamine Standard 2x1 ml  
Human serum. Concentration is given on the label, expressed of DMF (deoximorpholinofructose) and µmol/l of glycated albumin.

### STORAGE AND STABILITY OF REAGENTS

The Reagents are stable until the expiry date show on the label at 2-8°C.

Indications of deterioration: Presence of particulate material, turbidity. Absorbance of control at 530(±20) nm ≥0,065.

### SAMPLES

Don't use haemolized samples. Fructosamine in serum is stable 7 days at 2-8°C.

### REFERENCE VALUES

Serum: 1,9 – 2,9 mmol/l DMF, 217 – 331 µmol/l glycated albumin. Concentrations are slightly lower 5% in child.

Fructosamine reference values depends on albumin concentration. Plasmas give lower results than serum.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer thermostable at 37°C, able to read at 530(±20)nm.  
Pipettes for 50 µl and 1,0 ml.

### PRECAUTION

The kit is only for in vitro use.  
However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### REAGENT PREPARATION

Reagent is provided ready to use.  
Working standard: reconstitute Fructosamine Standard with 1,00 ml of distilled water. Mix gently and incubate 30 minutes before using.

The solution is stable for 15 days at 2-8°C if contaminations are prevented during their use, or for 45 days at -20 °C when frozen in aliquots.

### NOTES

The reagent may be used in several automated analyzers. Instructions for many of them are available on request.

### PROCEDURE

#### Assay conditions

Method: kinetic, FXT  
nitroblue tetrazolium NBT  
Wavelength: 530 (± 20) nm  
Temperature: 37°C  
Blank: distilled water

1. Bring the Reagent to room temperature (16-25)°C

2. Pipette into labelled test tubes:

	Standard	Test
Reagent, ml	1.0	1.0
Standard, µl	50	-
Sample, µl	-	50

2. Mix thoroughly and incubate at 37°C. Start the stopwatch.

3. Read the absorbance (A) of the Standard and the Test at 530 nm after exactly 10 minutes (A<sub>1</sub>) and 15 minutes (A<sub>2</sub>) against distilled water.

### CALCULATIONS

The fructosamine concentration in the sample is calculated using the following general formula:

$$\frac{(A_2 - A_1)_{\text{Sample}}}{(A_2 - A_1)_{\text{Standard}}} \times C_{\text{Standard}} = C_{\text{Sample}}$$

### QUALITY CONTROL

It is recommended to use the normal and patologic serum control.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### METROLOGICAL CHARACTERISTICS

Detection limit: 0,14 mmol/l DMF, 16 µmol/l glycated albumin.

Linearity limit: 7 mmol/l DMF, 800 µmol/l glycated albumin. For higher values dilute sample ½ with distilled water and repeat measurement.

Repeatability (within run):

Mean Concentration	CV	n
3,9 mmol/l = 446 µmol/l	2,7%	20
5,7 mmol/l = 651 µmol/l	2,5%	20

Reproducibility (run to run):

Mean Concentration	CV	n
3,9 mmol/l = 446 µmol/l	4,3%	25
5,7 mmol/l = 651 µmol/l	4,0%	25

\* CV - coefficient of variation n – number of determinations

Sensitivity: 21,2 mA x l/mmol DMF = 0,17 mA x l/µmol GA.

Interferences: Hemoglobin (10g/l), bilirubin (20 mg/dl) and lipemia (triglycerides 10 g/l) do not interfere. Other drugs and substances may interfere.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

Fructosamine is the generic name for plasma protein ketoamines formed by non-enzymatic attachment of glucose to amino groups of proteins (largely albumin). The measurement of fructosamine is useful to monitor the average concentration of blood glucose for an extended period of time (2-3 weeks) in individuals with diabetes mellitus. Because fructosamine determination monitors short-term glycemic changes different from glycated hemoglobin, it is recommended to be used in conjunction with glycated hemoglobin rather than instead of it.

### BIBLIOGRAPHY

1. Baker R John, Metcaff A Patricia, Johnson N Roger, Newman David and Rietz Peter. Use of protein based standards in automated colorimetric determinations of fructosamine in serum. Clin Chem 1985; 31: 1550-1554.
2. Tietz Textbook of Clinical Chemistry, 2<sup>nd</sup> edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
3. Van Diejen-Visser MP, Seynaeve C and Brombacher PJ. Influence of variations on serum fructosamine concentration. Clin Chem 1986; 32: 1610.
4. Hurst L Paul. Effect of anticoagulants on fructosamine determination. Clin Chem 1987; 33: 1947.
5. Young DS. Effects of drugs on clinical laboratory tests 3th ed. AACC Press, 1997.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed AACC Press, 1997.