DAC-SpectroMed s.r.l.

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Hemo-DAC.Lq

HEMOGLOBIN

DRABKIN METHOD WITH HEXACYANOFERRATE

For « in vitro » use only
 Store at 2-8°C

Cod 3048H2000 2000 ml
 Cod 3048H3000 3000 ml

PRINCIPLE

Ferrous ions of hemoglobin are oxidized to the ferric state by potassium ferricyanide to form hemoglobin (methemoglobin). Hemoglobin reacts with cyanide to form hemoglobincyanide (cyanmethemoglobin) medium forming a colored complex. The intensity of coloration, measured at 540(±20) nm, is proportional to hemoglobin concentration^{1,2}.

CONTENT AND COMPOSITION

Concentrated Reagent (x 25)
 Sodium hydrogen carbonate 0,30 mol/l
 Potassium hexacyanoferrate (III) 15,2 mmol/l
 Acetone cyanohydrine 0,14 mol/l

Toxic! Do not pipete by mouth!

Never allow reagent to come in contact with acid!

Hemoglobin Standard 2 ml
 Hemoglobin ≅ 120 g/l. Concentration is given on the label.

STORAGE AND STABILITY OF REAGENTS

Reagents at 2-8°C are stable until the expiry date shown on the label when stored tightly closed, contaminations are prevented during their use and are protected from light.

SAMPLES

Venous or capillary blood collected by standard procedures with heparin or EDTA as anticoagulants.
 Hemoglobin in blood is stable for 6 days at 2-8°C.

QUALITY CONTROL

It is recommended to use commercial whole-blood control materials to verify the performance of the measurement procedure.
 Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 540±20 nm. Timer. Dropper for 10 µl and 2,5 ml.

PRECAUTION

For in vitro diagnostics only.
 Handle all patients' samples as potentially dangerous and treat as infectious.
 Precautions established for work with caustic and toxic substances should be observed while using the reagents.

REAGENT PREPARATION

Working Reagent:
 Dilute the contents of the Concentrated Reagent bottle up to 1000 ml with distilled water. Mix thoroughly.
 Other volumes can be prepared by mixing in the proportion:
 1 ml Concentrated Reagent + 24 ml distilled water.
 Store the diluted reagent in a brown glass bottle.
 Stable for 6 months at 15-30°C. Do not freeze it.
 Hemoglobin Standard provided ready to use

PROCEDURE

Assay conditions
 Method: end point
 Wavelength : 540 (±20) nm
 Light path: 1 cm
 Temperature : 15-25°C
 Blank: against reagent

1. Pipette into labeled test tubes:

	Blank	Standard	Sample
Hemoglobin Standard, µl	-	10	-
Sample, µl	-	-	10
Working Reagent, ml	2,5	2,5	2,5

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

2. Mix thoroughly and incubate the tubes for 20 minutes at room temperature (16-25°C).
 3. Measure the absorbance (A) of the Standard and the Sample at 540(±10) nm against the Blank.

The color is stable for at least 2 hours.

CALCULATIONS

The hemoglobin concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{St}} = C_{\text{Sample}}$$

REFERENCE VALUES

Men³ 130-160 g/l
 Women³ 120-140 g/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

Detection limit: 20 g/l.
 Linearity limit: 200 g/l.

Repeatability (within run):

Mean Concentration	CV	n
100 g/l	2,3 %	20
150 g/l	1,1 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
100 g/l	3,1 %	25
150 g/l	2,1 %	25

* CV - coefficient of variation n – number of determinations

Interferences: Lipemia may falsely elevate the hemoglobin concentration because of turbidity. Bilirubin does not interfere. Other drugs and substances may interfere⁴. These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Hemoglobin is a red-pigmented, oxygen carrying protein of the red blood cells of vertebrates. It is a tetramer consisting of two pairs of unlike polypeptide chains, carrying an iron-containing porphyrin derivative called heme. Hemoglobin values are influenced by physiological variations as age, sex, physical exercise, posture, dehydration and altitude, and by pathological processes. Pathological levels are encountered in anemias and in polycythemias. Impaired production of red cells by the bone marrow, excessive blood loss and impaired delivery of red cells to the peripheral blood are the three main causes of anemia^{3,5}. Increased hemoglobin values are found in polycythemia vera, erythrocytosis, dehydration, newborns, acquired or congenital cyanosis, chronic heart and lung disease, renal cysts and a number of erythropoietin producing tumors^{3,5}. Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. International Committee for Standardization in Hematology. Recommendations for reference method for hemoglobinometry in human blood (ICSH Standard EP 6/2: 1977) and specifications for international haemoglobincyanide reference preparation (ICSH Standard EP 6/3: 1977). *J Clin Pathol* 1978; 31: 139-143.
2. Van Kampen EJ and Zijlstra WG. Standardization of hemoglobinometry: the hemoglobincyanide method. *Clin Chim Acta* 1961; 6:538-544.
3. Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
4. Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
5. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.