DAC-SpectroMed s.r.l.

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Cod 30561100 2x50 ml

PRINCIPLE

Transferrin-bound ferric ions in the sample are released by guanidine and reduced to ferrous by means of hydroxylamine. Ferrous ions react with ferrozine forming a colored complex. The intensity of coloration, measured at 560(±20) nm, is proportional to iron^{12,3}.

CONTENTS AND COMPOSITION

Reagent A	2x40 ml	pH 4.0
Guanidine chloride		1.0 mol/l
Acetate buffer		0.4 mol/l
Hydroxylamine		0.6 mol/l
Reagent B	2x10 ml	
Ferrozine		40 mmol/l
Ascorbic acid		200 mmol/l
Iron Standard	5 ml	

Concentration is given on the label. Aqueous primary standard.

NB: Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label when stored tightly closed.

Indications of deterioration:

Presence of particulate material, turbidity, absorbance of the blank over 0,050 at 560(±20) nm (1 cm cuvette).

SAMPLES

Serum or heparinized plasma collected by standard procedures.

Iron in serum or heparinized plasma is stable for 7 days at 2-8°C.

REFERENCE VALUES

Serum and plasma⁴:

Men: 65-175 µg/dl = 11.6 – 31.3 µmol/l.

Women: 50-170 µg/dl = 9.0-30.4 µmol/l.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 560±20 nm. Dropper for 200 µl and 1,0 ml.

PRECAUTION

The kit is only for in vitro use. All specimens must be considered potentially hazardous and handled as infectious.

REAGENT PREPARATION

Iron Standard is ready to use.

Working Reagent: 1 ml Reagent B + 4 ml Reagent A.

Mix thoroughly. Stable for 6 months at 2-8°C.

PROCEDURE

Method: end point
 Wavelength : 560 (±20) nm
 Light path: 1 cm
 Temperature : 15-30°C
 Blank: distilled water

1. Bring the Working Reagent to room temperature (15-30°C).
2. Pipette into labelled test tubes:

	Reagent Blank	Sample Blank	Sample	Standard
Distilled Water, µl	200	-	-	-
Sample, µl	-	200	200	-
Iron Standard, µl	-	-	-	200
Reagent A, ml	-	1,0	-	-
Working Reagent, ml	1,0	-	1,0	1,0

3. Mix thoroughly and let stand the tubes for 5 minutes at 16-25°C.

4. Read the Absorbance of the Sample Blanks (A_{SB}) at 560 nm against the distilled water.

5. Read the Absorbance of the Sample (A_{Sample}) and of the Standard (A_{St}) at 560 nm against the Reagent Blank.

CALCULATIONS

The iron concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}} - A_{\text{SB}}}{A_{\text{St}}} \times C_{\text{Standard}} \times \text{Sample Dilution Factor} = C_{\text{Sample}}$$

METROLOGICAL CHARACTERISTICS

Detection limit: 4 µg/dl = 0,71 µmol/l iron.

Linearity limit: 1000 µg/dl = 179 µmol/l iron.

For higher values dilute sample 1/2 with distilled water and repeat measurement.

Repeatability (within run):

Mean Concentration	CV	n
103 µg/dl = 18,4 µmol/l	2,2 %	20
305 µg/dl = 54,6 µmol/l	0,7 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
103 µg/dl = 18,4 µmol/l	2,9 %	25
305 µg/dl = 54,6 µmol/l	2,2 %	25

Sensitivity: 88 mA x dl/µg = 4,86 mA x l/µmmol.

Interferences: Bilirubin does not interfere. Do not use hemolyzed or lipemic sera. Other drugs and substances may interfere⁵. These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Iron is distributed in the body in a number of different compartments: hemoglobin, myoglobin, tissues (mainly in liver, spleen, and bone marrow). Only 0.1% of total body iron is present in plasma. Serum iron concentration is affected by many physiological or pathological conditions. Day-to-day variation is quite marked in healthy people. Iron deficiency and iron overload are the major disorders of iron metabolism. However, altered iron metabolism is also related to a number of other diseases.

Serum iron is increased in hemochromatosis, in acute iron poisoning, in active cirrhosis or acute hepatitis and as a result of increased transferrin levels^{3,5}. Serum iron concentration is decreased in many but not all patients with iron deficiency anemia and in chronic inflammatory disorders. Measurement of serum iron should not be used as a test for identification of iron deficiency^{3,5}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Garcic A. Highly simple determination of serum iron using chromazurol B. Clin Chim Acta 1979; 94: 115-119.
2. Paris M, Benoit MO, Rigat B, Progon J. A manual method for the direct determination of serum iron using new chromogen: Ann Biol Clin 1986; 44: 511-516.
3. Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
4. Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACCC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACCC Press, 1997