MONARCH™
IOL DELIVERY SYSTEM
CARTRIDGE

ENGLISH PRODUCT INFORMATION	2 - 6
DANSK PRODUKTINFORMATION	7 - 11
NEDERLANDS PRODUCTINFORMATIE	12 - 16
FRANCAIS INFORMATIONS SUR LE PRODUIT	17 - 21
DEUTSCH GEBRAUCHSANWEISUNG	22 - 26
ΕΛΛΗΝΙΚΑ ΠΛΗΡΟΦΟΡΙΕΣ ΓΙΑ ΤΟ ΠΡΟΪΟΝ	27 - 31
ITALIAN INFORMAZIONI SUL PRODOTTO	32 - 36
NORSK PRODUKTINFORMASJON	37 - 41
PORTUGUÊS FOLHETO INFORMATIVO	42 - 46
RUSSIAN ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ.....	47 - 51
ESPAÑOL INFORMACIÓN DEL PRODUCTO	52 - 56
SVENSKA PRODUKTINFORMATION	57 - 61
УКРАЇНСЬКА ІНСТРУКЦІЇ З ВИКОРИСТАННЯ	62 - 66

MONARCH™ IOL DELIVERY SYSTEM CARTRIDGE

DESCRIPTION

The MONARCH™ IOL Delivery System consists of two parts. An autoclavable, reusable, titanium handpiece (**Figure 1**) and a sterile, single-use cartridge (**Figure 2**) is used for implanting foldable intraocular lenses into the eye following removal of the natural crystalline lens. The system provides a controlled means to reliably place Alcon qualified intraocular lenses (IOLs) into the capsular bag when a qualified combination of handpiece, cartridge, Alcon foldable IOL, and ophthalmic viscosurgical device (OVD) are used.

The cartridge is loaded by inserting the IOL into the opening in the back of cartridge after an OVD is applied to the inner lumen of the cartridge. The loaded cartridge is installed into the handpiece and the IOL is delivered through the cartridge nozzle.


Figure 1:
Design of MONARCH™
Handpiece

Figure 2:
Design of MONARCH™
Cartridge

HOW SUPPLIED

The cartridge is supplied dry and individually sealed in a polyester/TYVEK[†] material pouch, which is terminally sterilized by Ethylene Oxide (EO).

INTENDED USE

The MONARCH™ IOL Delivery System, with reusable handpiece and single-use cartridge, is intended for use by a trained ophthalmic staff to implant Alcon qualified foldable IOLs into the eye following removal of the natural crystalline lens. The system provides a means to insert the IOL into the eye in a controlled manner.

INDICATION

The MONARCH™ IOL Delivery System is for implantation of qualified Alcon Foldable IOLs. No unqualified lenses should be used with the MONARCH™ IOL Delivery System.

QUALIFIED COMBINATIONS FOR IOL IMPLANTATION

The Alcon MONARCH™ cartridges are qualified for use with compatible Alcon MONARCH™ handpieces for the surgical implantation of Alcon qualified foldable IOLs. The handpiece-cartridge combinations are specified in the respective MONARCH™ handpiece Instructions for Use. Alcon Foldable IOLs are qualified for use with an Alcon qualified delivery system (handpiece and cartridge) and ophthalmic viscosurgical device (OVD) combination. **For the qualified combinations for implantation (handpiece, cartridge, IOL, and OVD), consult the Instructions for Use for the particular Alcon foldable IOL model being implanted, which can be found at www.ifu.alcon.com.** The use of an unqualified combination may cause damage to the IOL and potential complications during the implantation process (e.g., corneal endothelial cell damage).

INSTRUMENTATION

The Alcon foldable lenses should be loaded with standard holding forceps having rounded, non-serrated blades, or the MONARCH™ loading forceps 560.01, or equivalent.

CONTRAINDICATIONS

There are no known contraindications to the use of the MONARCH™ IOL Delivery System when used as recommended.

WARNINGS

Complications with the MONARCH™ IOL Delivery System used in conjunction with cataract surgery and IOL implantation include the following: physical injury, postoperative inflammation, prion infection, toxic reactions (i.e., toxic anterior segment syndrome), irritation, corneal endothelial cell damage, corneal decompensation, infection (e.g., endophthalmitis, systemic infection), retinal detachment, corneal edema, pigment dispersion, tissue damage (e.g., capsular tissue damage, iris damage, and zonular damage), decreased or blurred vision, unexpected postoperative refraction, visual disturbance, IOL decentration, tilt, or dislocation, and elevated IOP.

CARTRIDGE PRECAUTIONS

1. A high level of surgical skill is required for intraocular lens implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.
2. DO NOT store cartridges at temperatures over 30 °C (86 °F).
3. Contents are sterile unless package is open or damaged. DO NOT use the cartridge if the sterility has been compromised or if the sterile package has been unintentionally opened before use.
4. DO NOT resterilize the cartridge by any method.
5. DO NOT reuse the cartridge. The cartridge is for single use only. Reuse of this single-use device may result in serious injury, such as serious infection (e.g., endophthalmitis and systemic infection).
6. Potential risks from reuse or reprocessing include: a damaged cartridge, a damaged lens, an unexpected delivery outcome, and contamination of the cartridge.
7. Discard the cartridge after use.

ORDERING INFORMATION

Contact Alcon Sales Representative or Customer Service at 1-800-TO-ALCON in the United States with the correct Product Order Number.

SERIOUS INCIDENT REPORTING

Any serious incident that may reasonably be regarded as device related should be reported to Alcon Laboratories, Inc.:

By Phone:

International - Contact the local country office.

United States – Call 1 (800) 757-9780

By Email:

In Europe: qa.complaints@alcon.com

Website:

<http://www.alcon.com/contact-us/>

Each handpiece and cartridge is identified by a lot number which provides traceability, and this information should be given to Alcon.

NOTE: In Europe, serious incidents must also be reported to the competent authority for medical devices of the appropriate Member State.

EXPIRATION DATE

Sterility is guaranteed until the use-by date unless the primary sterilization package is damaged or opened. The use-by date is clearly indicated on the outer box label of the carton. Any cartridge held *after* the use-by date should NOT be used.

PERFORMANCE CHARACTERISTICS

An evaluation of postmarket surveillance, clinical investigations, and clinical literature data support that the MONARCH™ IOL Delivery System has an acceptable safety and performance profile when used to implant qualified Alcon foldable IOLs. Adverse events rates associated with using the MONARCH™ IOL Delivery System are low while clinical outcomes are favorable. Clinical outcomes were consistent with the intended use of the MONARCH™ IOL Delivery System.

CLINICAL FEATURES of the MONARCH™ IOL Delivery System

The MONARCH™ IOL Delivery System (handpiece and cartridge) is qualified with Alcon foldable IOLs and has the following features:

- The MONARCH™ IOL Delivery System allows for the implantation of qualified Alcon foldable IOLs into the eye following removal of the crystalline lens.
- The MONARCH™ IOL Delivery System provides a controlled means to reliably place the qualified Alcon foldable IOL into the capsular bag.

DIRECTIONS FOR USE

Handpiece Preparation

Step 1. Thoroughly clean and sterilize the handpiece prior to the first use.

NOTE: The handpiece must also be cleaned and then sterilized after each subsequent use. The handpiece must be sterilized in accordance with guidelines provided (refer to **Handpiece Cleaning and Handpiece Sterilization** sections below).

Step 2. Inspect the handpiece prior to each use to ensure that it is not damaged.

- 2.1. Push the knob forward until the threads contact the barrel.
- 2.2. Turn the knob clockwise to advance the plunger.
- 2.3. Ensure the plunger moves smoothly through the entire delivery process until the knob touches the housing.

IMPORTANT: The use-life of an un-damaged MONARCH II/III handpiece, when used per the IFU, is indefinite. In the event the handpiece, particularly the plunger tip, appears damaged, bent, or incapable of proper use, **DO NOT** use the handpiece and contact Alcon.

Step 3. Retract the plunger of the handpiece completely for cartridge loading.

Preparation of MONARCH™ Cartridge Prior to Loading and Folding of Alcon Foldable IOL

Step 1. Use the MONARCH™ cartridge at Operating Room temperatures between 18° C (64° F) and 23° C (73° F).


Step 2. Inspect the cartridge pouch carefully for tears, cuts, punctures, or other signs that the sterile package has been opened or damaged. This device is sterile until the primary package is opened. DO NOT use the cartridge if the sterility has been compromised or if the sterile package has been unintentionally opened before use. If sterility is compromised, use another cartridge.

Step 3. Open the pouch using standard sterile procedures and transfer the cartridge to the sterile field.

Step 4. Inspect the cartridge nozzle tip, the inner cartridge lumen, and the rest of the cartridge. If damage, particulates, or deformation is noted, please use another cartridge. Proceed to Step 5 with system components that pass the above checks.

Step 5. Only use an Alcon OVD qualified for use with the MONARCH™ IOL Delivery System that has been allowed to come to the operating room temperature. **Completely** fill the cartridge with OVD, as shown in **Figure 3**, immediately prior to loading and delivery of the lens. For a list of Alcon qualified OVDs, consult the Instructions for Use for particular Alcon foldable IOL model being implanted. DO NOT attempt to load the lens without adequate OVD in the device.

Figure 3: Filling of MONARCH™ cartridge with OVD


Lens Loading and Folding for Single-Piece Alcon Foldable IOL with Standard Forceps

NOTE: For lens loading and folding of the Single-Piece Alcon foldable IOL with the MONARCH™ 560.01 loading forceps, see instructions supplied with forceps.


Step 1. Using holding forceps, grasp the lens by the optic edge and gently place the lens anterior side up into the back of the OVD-filled cartridge as shown in **Figure 4**. The lens should be inserted until the optic is a little more than half-way inside the cartridge. Use the holding forceps to gently push down on the lens, verifying that the lens is on the bottom surface of the cartridge.

Figure 4: View of IOL insertion into OVD-filled MONARCH™ cartridge


Step 2. Using holding forceps, take the trailing haptic, and gently fold the haptic onto the anterior side of the optic as shown in **Figure 5**.

Figure 5: View of IOL haptic folding in OVD-filled MONARCH™ cartridge

Step 3. As shown in **Figure 6**, slowly grip or push the optic edge to position the lens as far into the cartridge as the forceps will permit, while ensuring the lens remains on the bottom surface of the cartridge and the trailing haptic remains on the optic.

Figure 6: Positioning of IOL in OVD-filled MONARCH™ cartridge

Alternate Loading and Folding Technique for Alcon Multi-Piece IOL into MONARCH™ II Cartridge

Step 1. Use holding forceps to grasp the lens by the optic edge and gently place the lens anterior side up into the back of the OVD-filled cartridge as shown in **Figure 7a**. The lens should be inserted until it is centered with or slightly past the outline etched into the top of the cartridge. Lenses with a 6.5 mm optic diameter may need to be pre-folded for easy entry into the back of the cartridge. Lenses may be pre-folded by gently applying pressure to the edge of the lens while holding the central optic with forceps as in **Figure 7b**. The trailing haptic will extend from the proximal end of the cartridge. Position the trailing haptic to the left of the haptic post as shown in the detailed view in **Figure 7a**. This will allow the plunger to go past the haptic during the initial advancement of the plunger into the cartridge.


Figure 7a: Insertion of IOL into MONARCH™ II Cartridge


Figure 7b: Pre-folding of IOLs with 6.5 mm optic

Step 2. Verify the lens is positioned on the bottom surface of the cartridge. The haptic should be maintained to the left of the post when the lens is loaded correctly. Accurate positioning of the lens will decrease the potential for optic and haptic damage.

Insertion of Cartridge into Handpiece and Lens Delivery

NOTE: For instructions on how to insert the cartridge into the MONARCH™ handpiece and deliver the lens, please refer to the Instructions for Use for the corresponding MONARCH™ handpiece.

NOTE: Refer to IOL DFU for dwell time limitations for a particular IOL model.

HANDPIECE CLEANING

! IMPORTANT

- Thoroughly clean the handpiece immediately before **first use** and **after each subsequent use**.
- **Manual cleaning** should be performed immediately after each surgical procedure, followed by cleaning using an **automated washing procedure**.
- Cleaning instructions can be located in the respective MONARCH™ handpiece Instructions for Use, which can be found at www.ifu.alcon.com.

HANDPIECE STERILIZATION

! IMPORTANT


- Thoroughly sterilize the handpiece immediately before **first use** and **after each subsequent use**.
- Sterilization instructions can be located in the respective MONARCH™ handpiece Instructions for Use, which can be found at www.ifu.alcon.com.

**STAND ALONE SYMBOLS FROM ISO 7000 USED ON LABELING
(ISO 7000/IEC 60417* Title: Graphical Symbols for Use on Equipment)**

Symbol	Reference Number from ISO 7000 / IEC 60417*	Symbol Title / Explanatory Text
	1051	Do not re-use
	2608	Do not resterilize
	2607	Use-by date
	2501	Sterilized using ethylene oxide
	2492	Batch code
	2493	Catalogue number
	0434A	Caution
	3082	Manufacturer
	0533	Upper limit of temperature
	1641	Consult instructions for use or consult electronic instructions for use
	2606	Do not use if package has been damaged
	2497	Date of Manufacture
	3079	Open Here
	3707	Single sterile barrier system
	6049*	Country of manufacture

*This symbol is from IEC 60417.

ABBREVIATIONS or OTHER SYMBOLS USED ON LABELING

Abbreviation or Symbol		Explanatory Text
		Medical device
IOL		Intraocular lens
		Cartridge
		Handpiece
OVD		Ophthalmic viscosurgical device
		Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician
		Authorized Representative in the European Community

Alcon Laboratories, Inc.
6201 South Freeway
Fort Worth, Texas 76134-2099 USA

Alcon Laboratories Belgium
Lichterveld 3
2870 Puurs-Sint-Amands
Belgium

U.S. Patents: www.alconpatents.com

† Trademarks are the property of their respective owners.

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© 2020 Alcon Inc.