## DAC-SpectroMed s.r.l.

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PT MD 11-15796482-001:2003

## Phosph-DAC.Lq

PHOSPHORUS  
 PHOSPHOMOLYBDATE UV  
 For « in vitro » use only  
 Store at 2-8°C

Cod 3064P200 2x100 ml

### PRINCIPLE

Inorganic phosphate reacts in acid medium with ammonium molybdate to form a phosphomolybdate complex with yellow color.

The intensity of the color measured at 340(±10) nm is proportional to the phospholipids concentration in the sample.

### CONTENTS AND COMPOSITION

Reagent	2x100 ml	pH 7,55
Ammonium molybdate		0,40 mmol/l
Sulphuric acid		210 mmol/l
Detergents		
Phosphorus Standard	5 ml	
Phosphorus aqueous primary standard.		

### STORAGE AND STABILITY OF REAGENTS

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protect from light and contaminations prevented during their use.

Indications of deterioration: Presence of particulate material, turbidity. Blank absorbance at 340(±10) nm ≥0,54.

### SAMPLES

Serum or plasma. Free of hemolysis. Serum or plasma should be removed from the clot as quickly as possible to avoid elevation of serum phosphorus from hydrolysis or leakage of phosphate present in erythrocytes. Stability: 7 days at 2-8°C.

Urine (24h): collect the specimen into a bottle containing 10 ml of 10% v/v hydrochloric acid (HCl) to avoid phosphate precipitations. Adjust to pH2. Dilute the sample 1/10 with distilled water. Mix. Multiply the result by 10 (dilution factor). Stability: 10 days at 2-8°C.

### REFERENCE VALUES

Serum or plasma.

Adults: 2,5 – 5,0 mg/dl = 0,80-1,61 mmol/l.

Children: 4,0 – 7,0 mg/dl = 1,29-2,26 mmol/l.

Urine: adults: 0,4 - 1,3 g/24-h.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer thermostatic at 37/30/25°C, able to read at 340(±10)nm.

Pipettes for 10 µl and 1,0 ml.

### PRECAUTION

The kit is only for in vitro use.

However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### REAGENT PREPARATION

Reagents are ready to use.

### QUALITY CONTROL

It is recommended to use the Sera N (code 2055S5) and Sera P (code 2057S5) to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### PROCEDURE

Assay conditions

Method: end point  
 Wavelength : 340 (±10) nm  
 Temperature : 25/30/37°C  
 Blank: against reagent

1. Pipette into labeled test tubes:

	Blank	Standard	Sample
Reagent	1,0 ml	1,0 ml	1,0 ml
Phosphorus Standard	-	10 µl	-
Sample	-	-	10 µl

2. Mix thoroughly and incubate for 5 min. at reaction temperature.

3. Measure the absorbance (A) of the Standard and the Sample at 340 nm against the Blank.

### CALCULATIONS

The phosphorus concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} \times C_d = C_{\text{Sample}}$$

Urine:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} \times V_{24h} = C_{24h}$$

mg/dl x 0,323 = mmol/l

### METROLOGICAL CHARACTERISTICS

Detection limit: 0,07 mg/l

Linearity limit: 15 mg/dl.

If the results obtained were greater than linearity limit, dilute the sample ½ with NaCl 9 g/l water and repeat measurement.

Repeatability (within run):

Mean Concentration	CV	n
3,44 mg/dl	0,63%	20
5,84 mg/dl	0,64%	20

Reproducibility (run to run):

Mean Concentration	CV	n
3,45 mg/dl	0,72%	20
5,83 mg/dl	0,68%	20

\* CV - coefficient of variation n – number of determinations

Sensitivity: 1 mg/dl = 0,053 A.

Trueness: Results obtained using 50 sample were the following:

Correlation coefficient: r = 0,9938

Regression equation: y = 0,9902x + 0,0749.

Interferences: hemolysis and other drugs and substances may interfere.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

Phosphorus is an essential mineral for tissue bone formation and is required by every cell in the body for normal function. Approximately 85% of body phosphorus is found in bone and in teeth. Low levels of phosphorus, can be caused by hypervitaminosis D, primary hyperparathyroidism, renal tubular disorders, antacids or malabsorption.

High levels of phosphorus can be caused by diet, bone metastases, liver disease, alcohol ingestion, diarrhea and vomiting.

### NOTES

1. Calibration with the provided aqueous standard may cause a matrix related bias, especially in some analyzers. In these cases, it is recommended to calibrate using a serum based standard (Multi St-DAC, cod. 2051M5).

2. This reagent may be used in several automated analyzers. Instructions for many of them are available on request.

3. Use clean disposable pipette tips for its dispensation.

4. The majority of detergents contain phosphorus, that's why it is recommended to wash the labware with diluted nitric acid, and then flush with water.

### BIBLIOGRAPHY

1. Ferrell E.C. phosphorus. Kaplan A et al. Clin Chem The C. V. Mosby Co St Louis, Toronto. Princeton 1984; 1072-1074 and 418.
2. Daly JA, et al. Clin Chem 1972; 18 (3): 263-265.
3. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
4. Young DS. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
5. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC Press, 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC Press, 1995.