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Ure - DAC.Lq

UREA


ENZYMATIC COLORIMETRIC METHOD -
UREASE/SALICYLATE

For «in vitro» use only
Store at 2-8°C

Cod 3092U200 200 ml
Cod 3092U500 500 ml
Cod 3092U1000 1000 ml

PRINCIPLE

Urea in the sample originates, by means of the coupled reactions described below, a colored complex that can be measured by spectrophotometry^{1,2,3}.

CONTENTS AND COMPOSITION

Reagent A	pH 6,9
Sodium salicylate	62 mmol/l
Sodium nitroprusside	3,4mmol/l
Phosphate buffer	20 mmol/l
Reagent B	
Urease	> 500 U/ml
Reagent C	
Sodium hypochlorite	7 mmol/l
Sodium hydroxide	150 mmol/l

Irritant! Irritating to eyes and skin.

Urea Standard. Aqueous primary standard

Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

Reagents and Urea Standard are stable at 2-8°C until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Indications of deterioration: Presence of particulate material, turbidity, absorbance of the blank over 0,100 at 600(±20) nm (1 cm cuvette).

SAMPLES

Serum, plasma or urine collected by standard procedures.

Urea in serum or plasma is stable for 7 days at 2-8°C.

Heparin is recommended as anticoagulant. Urea in urine is stable for 3 days at room temperature if microbial growth is prevented.

Dilute urine 1/50 with distilled water before measurement.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 600±20 nm. Thermostatic water bath at 37°C. Dropper for 10 µl and 1,0 ml.

PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

REAGENT PREPARATION

Reagent C and Standard are provided ready to use.

Mix carefully contents of Reagent B vial before using.

Working Reagent: can be prepared in the proportion: 1 ml Reagent B +24 ml Reagent A. Mix thoroughly. Stable for 2 weeks at 2-8°C.

PROCEDURE

Assay conditions

Method:	end point
Wavelength :	600 (± 10) nm
Light path:	1 cm
Temperature :	16-25°C/37°C
Blank:	against reagent

1. Bring the Reagents to room temperature (15-25°C).

2. Pipette into labeled test tubes:

	Blank	Standard	Sample
Urea Standard, µl	-	10	-
Sample, µl	-	-	10
Working Reagent, ml	1,0	1,0	1,0

3. Mix thoroughly and incubate the tubes for 10 minute at room temperature (15-25°C) or for 5 minutes at 37°C.

4. Pipette:

Reagent C, ml	1,0	1,0	1,0
---------------	-----	-----	-----

5. Mix thoroughly and incubate the tubes for 10 minute at room temperature (15-25°C) or for 5 minutes at 37°C.

6. Read the absorbance (A) of the Standard and the Sample at 600 nm against the Blank. The colour is stable for 2 hours.

CALCULATIONS

The urea concentration in the sample is calculated using the following general formula:

$$A_{\text{Sample}} \times C_{\text{Standard}} \times C_{\text{factor dilution Sample}} = C_{\text{Sample}}$$

If the Urea Standard provided has been used to calibrate :
Serum and Urine
plasma

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times 60 = \text{mg/dl urea} \quad \times 3000 = \text{mg/dl urea}$$

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times 10 = \text{mmol/l urea} \quad \times 500 = \text{mmol/l urea}$$

REFERENCE VALUES

Serum, plasma⁴:

Newborn: 1,8-6,4 mmol/l

Adult:

18-60 ani: 2,1-7,1 mmol/l

60-90 ani: 2,9-8,2 mmol/l

>90 ani: 3,6-11,1 mmol/l.

15-39 mg/dl urea = 7-18 mg/dl BUN = 2,5-6,5 mmol/l urea.

Concentrations in the neonatal period are lower, and in adults over 60 years of age are higher than in adults. Concentrations also tend to be slightly higher in males than in females.

Urine⁴:

26-43 g/24-h urea = 12-20 g/24 h BUN = 428-714 mmol/24-h urea.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

Detection limit:

1,3 mg/dl urea = 0,60 mg/dl BUN = 0,21 mmol/l urea.

Linearity limit: 300 mg/dl urea = 140 mg/dl BUN = 50 mmol/l urea.

For higher values dilute sample 1/5 with distilled water and repeat measurement.

Repeatability (within run):

Mean urea concentration	CV	n
26 mg/dl = 4,3 mmol/l	1,6 %	20
86 mg/dl = 14,2 mmol/l	0,8 %	20

Reproducibility (run to run):

Mean urea concentration	CV	n
26 mg/dl = 4,3 mmol/l	2,4 %	25
86 mg/dl = 14,2 mmol/l	1,3 %	25

* CV - coefficient of variation n - number of determinations

Sensitivity: 8,6 mA.dl/mg = 0,143 mA.l/mmol

Interferences: lipemia (triglycerides 10 g/l) and bilirubin (20 mg/dl) do not interfere. Hemolysis (hemoglobin 2 g/l) and elevated ammonia interfere. Other substances and drugs may interfere⁵.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Urea is synthesized in the liver as a by-product of the deamination of amino acids. Its elimination in the urine represents the major route for nitrogen excretion. Elevated urea concentration in plasma is found as a result of a high-protein diet, increased protein catabolism, after a gastrointestinal hemorrhage, mild dehydration, shock and heart failure or treatment with glucocorticoids (pre-renal uremia)^{4,6}.

Post-renal uremia is caused by conditions that obstruct urine outflow: nephrolithiasis, tumor or prostatic hypertrophy. The usefulness of urea as an indicator of renal function is limited by the variability of its plasma concentration as a result of no renal factors^{4,6}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Chaney AL, Marbach EP. Modified reagents for determination of urea and ammonia. *Clin Chem* 1962; 8:130-132.

2. Searcy RL, Reardon JE, Foreman JA. A new photometric method for serum urea nitrogen determination. *Amer J Med Technol* 1967; 33:15-20.

3. Tabacco A, Meattini F, Moda E, Tarli P. Simplified enzymic/colorimetric serum urea nitrogen determination. *Clin Chem* 1979; 25:336-337.

4. Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.

5. Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.

6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997